

Product Name :
Vibrating Machine

Product Code :
NLAB-ENGINEERINGLB89001



Description :

Vibrating Machine

Technical Specification :

Since our establishment, we are engaged in offering optimal quality Vibrating Machine that conforms.

Especially designed to perfection, it consists of frame mounted on four coiled spring that helps to carry a cube mould and an eccentric revolving shaft mounted in ball and roller bearing. Ideal for the operation on 220V, 50Hz, single phase power supply, this has vibration frequency of 12000± 400 R.P.M. per minute.

Rugged construction, negligible maintenance, high operational fluency, excellent performance and longer working life are some of the salient features of the Vibrating Machine.

Accessories:-

Flat Belts

Set of Springs.

KEW Cube Mould 70.6MM (I.S.I. Marked)

Lubricant - 100cc

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003