

**Product Name :**  
Vibrating Machine

**Product Code :**  
NLAB-ENGINEERINGLB89001



**Description :**

Vibrating Machine

**Technical Specification :**

Since our establishment, we are engaged in offering optimal quality Vibrating Machine that conforms.

Especially designed to perfection, it consists of frame mounted on four coiled spring that helps to carry a cube mould and an eccentric revolving shaft mounted in ball and roller bearing. Ideal for the operation on 220V, 5Hz, single phase power supply, this has vibration frequency of  $12000 \pm 400$  R.P.M. per minute.

Rugged construction, negligible maintenance, high operational fluency, excellent performance and longer working life are some of the salient features of the Vibrating Machine.

Accessories:-

Flat Belts

Set of Springs.

KEW Cube Mould 70.6MM (I.S.I. Marked)

Lubricant - 100cc

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003