

**Product Name :**  
Video Microscope System 1500x

**Product Code :**  
NLAB-EDUSCHOOLEQ173005

**Description :**

Video Microscope System 1500x

**Technical Specification :**

Stand: Large base boom stand with 20mm mounting post. Base is 15" x 15.5". Height of vertical post is 17.5".

Camera Mount: Includes c-mount adapter to mount camera directly to the lens system.

Focusing: Includes coarse and fine focusing.

Working Distance: 7.9mm

Magnification: When using 19" monitor on-screen magnification is 1507x. Larger monitor = more magnification, smaller = less.

Camera: 3.2 Megapixel digital color camera with 1/2" CMOS sensor. Captures images 2048x1536. Captures motion video at 45fps at 680x510, 32fps @ 1048x768, 11.5fps @ 2048x1536. Includes basic capture and measurement software.

Illumination: Single arm light guide with high intensity LED illuminator, shining through objective.

Objective Lens: Your choice of the following lens:

Plan Achromat LWD 50x Objective, NA 0.55, WD 7.9mm, resolution 0.61µm

Semi Brightfield Metallurgical 50x Objective, NA 0.80, WD 1.2mm, resolution 0.42µm.

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com