

Product Name :
Visual Measurement Microscope

Product Code :
NLAB-EDUSCHOOLEQ173012

Description :

Visual Measurement Microscope

Technical Specification :

Eyepieces: Includes 10x eyepieces with 25mm reticle installed. Binocular head, eyetubes are inclined 35° and head rotates 360°.

Focusing: Rack and pinion focusing controls, 120mm focusing stroke.

Stereo Microscope Body : Zoom ratio of 1:6.7. Greenough system designed for harsh environments of production and testing. Body has detent for duplicating exact measurements at different magnification settings.

Total Magnification Range options: 7.5 - 50x.

Illumination: Your choice of LED ring light or LED dual pipe light.

Stand: Track stand includes two stage clips. Stand is ergonomic, with low-profile base. Includes black/white stage plate.

Reticle: Measuring reticle comes installed in eyepiece. Your choice of 10mm / 100 division reticle or 0.4" / 100 division reticle.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003