

Product Name :
Vital Sign Monitors

Product Code :
NLABMEDICAL-M3360008

Description :

Vital Sign Monitors Superior India, Vital Sign Monitors Superior manufacturer India, suppliers, Medical Vital Sign Monitors India, Hospital Vital Sign Monitors India

Technical Specification :

7"Color TFT Screen, waveform up to 8 channels Light and portable with build-in rechargeable battery 72 hour storage and review of trend grams and tables.

40 second review of holographic waveform RMS.

All-round monitor of adult, pediatric and neonatal electrosurgical unit and defibrillation protected.

Standard parameters: ECG, NIBP, SpO2, Respiration, Temperature, Pulse rate.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003