

| | |
|---------------------------------------|--|
| Product Name : Voltaic Cell | Product Code : NLAB-EDUSCHOOLEQ1730203 |
|---------------------------------------|--|

| |
|--|
| |
|--|

Description :

Voltaic Cell

Technical Specification :

The includes one 4" copper and one 4" zinc electrode for use in classroom experiments.

Ten ounce plastic cup to contain electrolyte solution is fitted with metal clamps to hold the electrodes in place.

The voltaic Cell with 2 Electrodes Voltaic cell demonstrates the generation of electricity through oxidation and reduction reactions for students in grades 9-12.

It comes with instructions.

The set includes one 4" copper and one 4" zinc electrode for use in classroom experiments.

Dimension:-

Width x Height: 99.06 mm x 96.52 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003