

**Product Name :**  
Voltmeter for Lab

**Product Code :**  
NLAB-EDUSCHOOLEQ94001



**Description :**

Voltmeter

**Technical Specification :**

- For demonstrating decomposition of water.
- This design is preferred for work of best precision.
- Platinum electrodes are used electrolysis of water.
- The electrodes are mounted with rubber stoppers for insertion in the open ends of the voltmeter limbs.
- The unit has two connected limbs, integral with reservoir tube and funnel shaped bulb, with a stop clock at the top of each limb.
- With suitable metallic stand, helps to prevent accidental breakage.
- Glass parts are made from corning glass

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003