

Product Name :
Ward Locker 18 Cabinet

Product Code :
NLABMEDICAL-M2570013



Description :

Ward Locker 18 Cabinet

Technical Specification :

Specification :-

Ward Locker is made of MS Sheet.

It has 18 Lockers with Individual Locks.

Mild Steel & Pre Treated Epoxy Powder Coated

We are Manufacturer, Exporters and Suppliers of Ward Locker 18 Cabinet at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003