

<b>Product Name :</b> Wash Basin Stand	<b>Product Code :</b> NLABMEDICAL-M2730013
<b>Description :</b> Wash Basin Stand	
<b>Technical Specification :</b> Height 840mm two Stainless steel bowl of 35 mm. Dia. Tubular ms steel frame work mounted on three swivel castors 50mm. Dia. Finish in epoxy powder coated	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003