

Product Name :
Water Bath

Product Code :
NLAB-SCIENTIFICEQ34002



Description :

Water Bath

Technical Specification :

Thermostatic Water Baths

Use for precision constant temperature and auxiliary.

Characteristics:

Inner tank: 1 hole, 2 holes, 2 rows & 4 four holes, 2 rows & 6 holes are stainless steel stretching inner chamber.

1 row & 4 holes, 1 row & 6 holes, 2 rows & 8 holes are stainless steel inner chamber.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003