

Product Name : Water Cooler	Product Code : NLAB-VOCATIONALTVET27097
Description :	
Water Cooler	
Technical Specification :	
Number of Taps : 2	
Dimensions : 24x22x50	
Cooling Capacity : 40 L/Hr	
Weight : 70 kg	
Material: Stainless Steel	
Storage Capacity : Greater Than 120 L	
Power Consumption (In Watts) : 800 - 1000 W	

Material Grade: SS 302

Features:

Reliable

Durable

Optimum quality

Fine designs

Sturdy

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

