

**Product Name :**  
Water Demineralizer, Type DI-425

**Product Code :**  
NLAB-EDUCATIONTVET13008



**Description :**

Water Demineralizer, Type DI-425

**Technical Specification :**

No pressure

Flow: 50 l / h

Capacity: 425 liters (@ 10 ° total salt content)

Water purity: 0.1 to 20  $\mu\text{S} / \text{cm}$

Connection: R 3/4 "

Supplied with analog conductivity meter (0.1 to 50  $\mu\text{S} / \text{cm}$ ), wall support and set of pipes.

Replacement cartridge for DI-425

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)