

Product Name :
Water Distiller Type GFL 2001/2

Product Code :
NLAB-EDUCATIONTVET13009



Description :

Water Distiller Type GFL 2001/2

Technical Specification :

Stainless steel heating body

Security lack of water

Capacity: 2 l / h

Conductivity: approx. 2.3 $\mu\text{S} / \text{cm} @ 25^\circ \text{C}$

Water consumption: approx. 20 l / h

Power: 2 kW

Power supply: 220V - 50Hz

Dimensions (WxDxH): 280 x 250 x 490 mm

Weight: 7.5 kg

Connection kit, for GFL 2001

Includes: inlet and outlet hoses (1.5m) as well as hose clamps

Plastic container, 5 liters

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003