

**Product Name :**  
Water Distiller, Type Gfl 2004

**Product Code :**  
NLAB-EDUCATIONTVET12011



**Description :**

Water Distiller, Type Gfl 2004

**Technical Specification :**

With 8 liter tank

Stainless steel tank and heater

Security lack of water

Capacity: 4 l / h

Conductivity: approx. 2.3  $\mu\text{S} / \text{cm}$  @ 25 ° C

Water consumption: approx. 48 l / h

Power: 3 kW

Power supply: 220V - 50Hz

Dimensions (WxDxH): 620 x 330 x 460 mm

Connection kit, for GFL 2002 to 2012 (includes: inlet and outlet hoses (1.5m) as well as cable ties

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003