

Product Name :
Water Distiller, Type Gfl 2004

Product Code :
NLAB-EDUCATIONTVET12011



Description :

Water Distiller, Type Gfl 2004

Technical Specification :

With 8 liter tank

Stainless steel tank and heater

Security lack of water

Capacity: 4 l / h

Conductivity: approx. 2.3 μ S / cm @ 25 ° C

Water consumption: approx. 48 l / h

Power: 3 kW

Power supply: 220V - 50Hz

Dimensions (WxDxH): 620 x 330 x 460 mm

Connection kit, for GFL 2002 to 2012 (includes: inlet and outlet hoses (1.5m) as well as cable ties

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003