

Product Name :
Water Multimeter

Product Code :
NLAB-SCIENTIFICIND13003



Description :

Water Multimeter

Technical Specification :

Multimeter for DO, Conductivity, and pH Multi Dual Input Portable Meter for pH & Conductivity or pH & DO
Standard and rugged (field) electrodes available.

500 measurement point storage memory,

Simple data transfer to printer or computer,

INTELLICAL electrodes that automatically remind user when calibration is necessary,

Resolution (selectable): 0.1/0.01/0.001

Simple calibration with screw on calibration sleeves,

Measuring range pH: 0 – 14,

Measuring range Conductivity: 0.01 μ S/cm – 400 mS/cm,

Measuring range Oxygen (LDO): 0.00 - 20.0 mg/l; 0 - 200% For measuring pH, dissolved oxygen conductivity and temperature of water samples while in the field (lake, streams, ponds, etc)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003