

Product Name :
Water Permeability Test Kit

Product Code :
ConcreteLab0004

Description :

Water Permeability Test Kit

Technical Specification :

Pressure chamber unit with 0-1.5 bar gauge & built-in micrometer gauge & valves.

Wrench for pressure lid.

Extra 0-6.0 bar gauge.

Water filling cup.

Adjustable clamping pliers.

Set of anchoring tools with drift & hammer.

Wrenches: 14 and 17 mm.

Sealant tape.

Bottles with boiled water, 3 bottles.

Gaskets, 10 mm thick, 4 pcs.

Gaskets, 15 mm thick, 4 pcs.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773