

Product Name :
Water Purification System

Product Code :
NLAB-SCIENTIFICIND13008



Description :

Water Purification System

Technical Specification :

Millipore Direct purifier System for ultrapure (Type I) water for analytical and biological laboratories from tap water.

Capable of removing ionic contaminations and Total Organic Carbon content (TOC) in water.

240V/50 Hz power supply

For purified water for reagent preparation.

High Flow/Capacity (4-10 Liters/Min) versions, for operation on either Ordinary Tap Water or Pretreated Water.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003