

Product Name :
Water Retention Apparatus

Product Code :
NLAB-CEMETLAB130026



Description :

Water Retention Apparatus

Technical Specification :

Two models available:

one fitted with aspirator pump and analog vacuum gauge with regulator

Water retention apparatus including an aspirator pump, vacuum regulator, vacuum gauge, three way stopcock, flask, rubber gasket, glass funnel, perforated brass dish and stand.

Weight: 8 kg (approx.)

the other with a portable vacuum pump and digital vacuum gauge with regulator

Water retention apparatus including vacuum pump with ultimate vacuum 0.1mbar, digital vacuum regulator complete with vacuum gauge 1mbar resolution, three way stopcock, flask, rubber gasket, glass funnel, perforated brass dish and stand.

230 V/50-60 Hz/1ph

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003