

Product Name :
Water Still Table Model

Product Code :
NLAB-LABORATORYEM31002



Description :

Water Still Table Model

Technical Specification :

Water Still Table Model

Double wall construction having both made out of highly polished stainless steel sheet.

On 220 V AC 50 Hz single phase.

Output Approx.

4/5 Liters per hour

25/26 Liters per hour

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003