


Product Name : Water Testing Lab Equipments	Product Code : NLAB-SCIENTIFICEQ32001
 A digital water testing instrument, likely a pH meter or conductivity meter, is shown. It has a black body with a small LCD screen displaying '9.39' and '23.0'. A black probe is held by a clamp and is submerged in a beaker containing a yellow liquid. The beaker is labeled '150 ml'.	
Description :	
Water Testing Equipments, Water Testing Lab Equipments	
Technical Specification :	
Water Testing Equipments - Hanna COD Analyzer, PH Meter, Distillation Unit and Dissolved Oxygen Meters offered by SGM Lab Solutions	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003