

Product Name : Water Vapour Permeability	Product Code : NLAB-TESTINGLABM45003
Description :	
Water Vapour Permeability	
Technical Specification :	
<p>Water Vapour Permeability Tester made of twelve brass bottles. It measures water vapour permeability of leather. The bottles have nut of circular opening. The bottles are mounted on the bottles holders which rotate at specified revolution per minute with the help of motor pulley & belts. The axis of the bottles is parallel to the axle of the bottles holders. A fan mounted in front of the mouths of the bottles and consisting of three blades in planes that is inclined at 120° to one another. The planes of the blades pass through the prolongation of the axle of the wheel. A pre set Digital timer is mounted for setting total no. of hours of instrument movement. All the components are chrome/ Zinc Plating & Painting for corrosion resistance long life.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003