

Product Name :
Welders Apron

Product Code :
NLAB-VOCATIONALTVET13033



Description :

Welders Apron

Technical Specification :

Lining: Unlined.

Product Style: Leather.

Size: One size fits most.

Color: Brown.

Material Type: Split cowhide.

Superior protection.

Heavy duty stitching.

Reinforced snaps

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003