

Product Name :
Welders Apron

Product Code :
SWEL-IN-0021NLAB-
MACHINELABWK34023



Description :

Welders Apron

Technical Specification :

Brown

Size: 107 x 60 cm

With flame retardant

Standard Welding Equipment Laboratory Set Welders Apron Manufacturer, Supplier and Welders Apron Exporter in India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

