

Product Name :
Welders Leather Gauntlet Gloves

Product Code :
SWEL-IN-0021NLAB-
MACHINELABWK34022



Description :

Welders Leather Gauntlet Gloves

Technical Specification :

Strong, Soft and supple leather

Soft lining

Sizes: Medium and large; Overall length 340mm

Standard Welding Equipment Laboratory Set Welders Leather Gauntlet Gloves
Manufacturer, Supplier and Welders Leather Gauntlet Gloves Exporter in
India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

