

Product Name :
Wet Sieve Shaker, Hand Operated

Product Code :
NLAB-TESTINGLABM109074



Description :

Wet Sieve Shaker, Hand Operated

Technical Specification :

Wet Sieve Shaker, Hand Operated

It consists of a motor with reduction gear box mounted at the top of the reservoir. The sieve holder, capable of accommodating seven sieves of 150 mm or 200 mm diameter is hung from an eccentric which imparts up and down motion to the nest of sieves held in the sieve holder and immersed in water. Supplied without sieves, and lid & receiver.

Wet Sieve Shaker, Motorized

Same as above, except that it can accommodate four sets of seven sieves each simultaneously for wet sieving. Supplied without sieves, and lid & Receiver. Suitable for operation on 220 V, 50 Hz, single phase supply.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003