

Product Name :
Wide Wheel Abrasion Testing Machine

Product Code :
NLAB-AGRETLAB160008



Description :

Wide Wheel Abrasion Testing Machine

Technical Specification :

Abrasion Testing Machine is designed for determining the resistance to abrasion/wear of natural stones and concrete products used for paving. The abrasion wheel is 70 mm thick and rotates at a speed of 75 rpm. The machine is equipped with a digital counter which stops the machine at the end of a preset number of revolutions.

Narrow Wheel

Abrasion Testing Machine is designed for determining the resistance to abrasion of unglazed tiles, grouts used for tiles and clay pavers. The abrasion wheel is 10 mm thick. Abrasive corundum sand, calibration marble (boloneisse) and mould (for EN 12808-2) should be ordered separately.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003