

Product Name :
Wire Gauge

Product Code :
TN391



Description :

Wire Gauge

Technical Specification :

Wire Gauge

Standard wire gauge measure either SWG/AWG standard

Round AWG SWG Wire Thickness Measuring Tester Ruler Gauge Diameter Steel Tool 1 Piece

Material: Stainless Steel

Color: Silver.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003