

Product Name :
Wood Turning Lathe Bench Model

Product Code
NLAB-MACHIN




Description :

Wood Turning Lathe Bench Model

Technical Specification :

Wood Turning Lathe Bench Model .

switch spur and
cup face plate with screw centre.



Length of lathe bed

3000mm

Spindle bearings

31020, 31020, 31020

660 RPM

Distance between centres

Spindle taper

260mm

Total weight

HP recommended

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003