

Product Name :
Wooden Tripod Heavy

Product Code :
NLAB-ENGINEERINGLB91002



Description :

Wooden Tripod Heavy

Technical Specification :

We provide our clients with wooden telescopic tripods that have been manufactured using the finest quality of wood.

This Wooden Tripod is used for providing stable support to surveying instruments such as theodolites, total stations, levels, and transits.

The wooden telescopic tripods have gained wide recognition for their durability and strength.

Our team of expert quality professionals carries out numerous stringent tests and checks on each product, so as to maintain international standards.

The wooden telescopic tripods can be customized keeping in mind the demands of our clients.

Features:

Provides stable support to surveying instruments such as theodolites, total stations, levels and transits

Manufactured using the finest quality of wood.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003