

**Product Name :**  
Wooden Tripod

**Product Code :**  
NLAB-ENGINEERINGLB91001



**Description :**

Wooden Tripod

**Technical Specification :**

We are recognized as the foremost manufacturers and suppliers of Wooden Telescopic Tripods which are utilized to support surveying instruments like theodolites, levels, transits and total stations. Our complete product range is manufactured from high grade raw materials and is available in a power coated finish.

Appreciated for features like high strength, fine finish and robust construction.

These Wooden Telescopic Tripods are widely demanded in the market.

In addition to this, offered tripods are available in various dimensions and other related specifications in order to fulfill the diversified requirements of our esteemed clients.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003