

Product Name :
Work Position Harness Nep

Product Code :
NLAB-SAFETYLABP20030



Description :

Work Position Harness Nep

Technical Specification :

Special Features: H-shaped thigh strap holds the positioning belt tight & fix makes it suitable for self rescue and ascending.

Anchor point: Work positioning waist belt with Two lateral.

D-rings and one ventral positioning D-ring.

Easy work positioning waist belt and sit harness

Ring provided for tool holder

Thermoform comfort back and thigh support

Materials: Fray-proof dope-dyed polyester webbing

Two lateral D-rings of forged alloy steel

Stitching: High-strength polyester stitching

Weight: 1.370 kgs approx.

Adjustability: Size adjustment straps at waist and thighs

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003