

**Product Name :**  
Work Position Harness NP 2001

**Product Code :**  
NLAB-SAFETYLABP20051



**Description :**

Work Position Harness NP 2001

**Technical Specification :**

Materials: Fray-proof dope-dyed polyester webbing Two lateral D-rings of forged alloy steel.

Anchor Point: Work positioning waist belt with Two lateral D-rings.

Special Features: Ring provided for tool holder Thermoform comfort back support.

Stitching: High-strength polyester stitching.

Adjustability: Size adjustment straps at waist.

Weight: 0.640 kgs approx.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

