

**Product Name :**  
Workbench 200Kg

**Product Code :**  
WC415



**Description :**

Workbench 200Kg

**Technical Specification :**

Workbench  
Commercial-grade steel frame with padded leveling feet.  
Easy to assemble  
Scratch-resistant  
Dimensions (Length x Width x Height) 1200 x 700mm x700mm.  
Without under storage  
H-frame white Laminated play wooden  
Max load 200Kg  
2 halls for computer cable passage

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003