

Product Name :
X Ray Machine-300 500 MA

Product Code :
NLABMEDICAL-M3350008

Description :

X Ray Machine-300 500 MA

Technical Specification :

Avis-300/500 is the new revolutionary range of X-ray

Its Image Quality and Contrast Help easy diagnosis.

It has been developed by keeping in mind all the problems faced by X-ray Developers.

It is easy to operate.

Avis 300/500 generator is design to fulfill the requirements of all Radiography Fluoroscopy

Spot Filming and special Procedures with manual or motorized table.

The generator is suitable for use with Image intensifier TV System.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003