

Product Name :
Yarn And Carpet Testing Lab Equipment

Product Code :
NLAB-TESTINGLABM48001



Description :

Yarn And Carpet Testing Lab Equipment

Technical Specification :

Readability 001g

Base size 120mm

Display Liquid Crystal Display (LCD) with Backlight

Operation AC Adapter (Included)

Communication RS232

Construction Metal base, ABS top housing, stainless steel pan,
glass draft shield with sliding top door, replaceable in-user cover

Unit tex, den, Nm, g, mg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003