

Product Name :
Yarn Evenness Tester

Product Code :
NLAB-TESTINGLABM48003

Description :

Yarn Evenness Tester

Technical Specification :

Yarn Evenness Tester is a kind of Yarn Testing Equipment which is used to measure the yarn appearance on boards. Manual and motorized models are available. The system consists of a yarn wrap board and an arrangement to wrap the yarn over it at equal spacing under a desired tension. The board is rotated manually to wrap yarn over it; the yarn is guided on the board through a yarn guide which moves parallel to the pivoting axis of the board with the help of a screw giving uniformly spaced wraps. The guide is positively linked to the rotation of the board with the help of a chain and sprocket drive. A set of change gears gives five different spacings to the yarn for wrapping yarn of different counts. A spring-loaded adjustable tensioning arrangement provides the desired tension to the yarn. The equipment is built on a rigid castled base. It is finished in dark gray painting and bright chrome/zinc plating to give it a corrosion-resistant finish. Two numbers of rectangular acrylic boards, having one face white and the other face black, to enable evaluation of both light and dark colored yarns are supplied as standard accessories. Additional boards

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003