

Product Name :
Youngs Modulus Apparatus

Product Code :
NLAB-MECHANICALAB240031



Description :

Youngs Modulus Apparatus

Technical Specification :

For determining the relationship between a load applied to a wire and the resulting extension. Consisting of an engraved brass scale with a vernier reading to 0.1mm. Complete with tension weight. Consists of a scale plate carrying a 0 to 30mm scale and a movable vernier readable to 0.1mm precision. Both have bars with clamping screws for the wires and hooks for the tension weight and loading masses. A ceiling clamp is provided for the upper ends of the wires, together with a pair of large wood screws for attaching it to a convenient overhead beam, door frame etc. A tension weight, mass 1.3kg approx..

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003